

SUPERVISED WIRELESS

EM-900 MONITOR STATION

Monitors up to **40+*** patients, people, locations and events

Wireless Long Range 900 MHz- up to 1000 ft. indoors- through walls, floors and ceilings!



*40+ Stations Monitored

The EM-900 sounds the Alert and LCD displays information such as:

123	J DOE	PX
-----	-------	----

Location Name Event

You may enter any information up to 15 characters

It's Portable! Can be carried with you! Alerts you to patient calls and emergency events while you're away from Nurse station or desk!

Now...
Advanced Wireless Technology provides **Fall Prevention, Nurse Call, Patient monitoring and Staff and Security notification** when needed.



Fall Prevention
Patient attempts to stand or walk


The **NEW EM-900 Monitor Station** is portable and monitors up to **40** patients at **Long Range!** The LCD Displays the alpha-numeric information that you have entered — received from wireless pendant and universal transmitter sensors. When sensors are activated, **alarm sounds** at Monitor Station and **LCD displays** information. Special features include **"Call Waiting"** for multiple alarms and data storage.



Nurse Call
Wireless & Convenient



"Walk Out"
Patients



Patients "At Risk"
to fall



Call for Help
Medical Staff and Security



Restricted
Areas

Now...

You have complete freedom to roam and monitor ALL of your patients!

How It Works:

Carry Monitor Station in hand, pocket, waist or shoulder strap!



EM-900 shown in re-charge cradle. Ni-Mh battery lasts up to 8 hours.

1. The EM-900 Monitor Station is **fast and easy** to set up and use. Enter information for each transmitter, ex: *Event - Name - Location*. You can enter any information (up to 15 characters) that will be displayed on LCD screen. It will *Beep Alert and Display* entered information whenever a *transmitter is activated*. **"Call Waiting"** will retain up to **10 multiple incoming calls**. Data Log retains calls for evaluation. Aux. out for Automatic Voice Dialers calls *cellphones, pagers and telephones*.

**Unlimited number of transmitters may be assigned to each of the 40 stations, however "missing transmitter" signal may be lost.*

PENDANT TRANSMITTERS



Pendant transmitters come complete with all accessories: Attachment cord, wrist band, necklace, surface mounting and belt clips.

Pendant Transmitters are attached to a cord, on wrist, necklace, beltclip or mounted. A **simple press of the button** sends the signal to EM-900 to sound alert and display LCD entered information. Pendants are supervised by EM-900 and **Automatically send alert** and displays whenever Pendant is missing such as a **"Walk Out"** patient. **Man Down** Pendant is Automatic and signals when patient is horizontal (Fall). All EM-900 transmitters **automatically** send signal to alert **low transmitter battery**.

UNIVERSAL TRANSMITTERS



Attached to bed and chair posts with Velcro straps or quick mount to surfaces.

2. Universal transmitters connect to sensor pads for Bed, Chair & Floor **Fall Prevention**. They also may be used with most types of sensors such as door switches, floor mats, temperature, IV sensors etc. When **activated** by sensor, Universal transmitters sends wireless signal to **sound local alert and display** LCD entered information. Instant and delay selection allow Fall Prevention patients **movement in bed or chair**. Transmitters turned off while attending patient. Missing transmitter and low battery signals **automatically sent**.

EM-900 transmitter range is up to 1,000 feet indoors, almost 1 mile open air & up to 5 miles with Repeater Unit!

SENSOR PADS



3. Sensor Pads activate Universal Transmitter when patient **attempts exit** from beds, wheelchairs or any place where **Fall Prevention** is needed. Pads also are used to detect and signal patient presence when walking in areas. *Many models are available*; Sensor strips for edge of bed and rails to large area non-slip corrugated models. Leads snap lock into transmitters **prevent patient removal**. Lifetime Sensor Pads are guaranteed to last a lifetime, eliminating replacement and expensive costs.

See Emergency Caller Product index Price list for all accessories for Products and Systems

Order Information

Model	Description
EM-900	Monitor Station complete with Monitor, Charging Cradle, AC Adapter (Order transmitters and sensors separately)
LP-900	Pendant Transmitter complete with Bed-Chair Cord, Wrist Band, Belt Clip, Necklace, Mount Clip, Batteries
LT-900	Universal Transmitter complete with Mounting Clip, Velcro Strap and Batteries
AT-900	Man Down Transmitter complete with Safety Pin, Necklace and Batteries
Sensor Pads	Please see Emergency Caller Lifetime Sensor Pad Catalog page for extensive listing of models
AVD-45be	Automatic Voice Dialer
PP-1	24 hour battery back up pack
LR-900	Repeater Transmitter Unit - Re-Transmits All 900MHz Transmitter Signals Up to 5 Miles

Emergency Caller Products



Emergency Caller Products (T) 800-227-2474 (F) 800-762-7613

www.emergencycaller.com sales@emergencycaller.com

U.S.A. and Foreign Patents Issued and Pending

MADE IN U.S.A.

